



Percutaneous Endoscopic Gastrostomy (PEG) Tube / Kit

Brand new design

- According to the feedback from market, Fushan redesigned the PEG kit totally from components to package. The new kit is more convenient to handle.

- Both economic pack and luxury pack PEG kit are available for choice.



commitment

- provides products with excellent quality and reliable service to meet the requirement of market.

Excellent Quality

- insists in manufacturing products with the highest standards. The Features of PEG kit include:

- Medical Grade Silicone Construction
- External Retention Ring
- Universal Feeding Port Adapters
- Collapsible Internal Retention Bumper
- Radiopaque Stripe and Bumper
- Tubing Clamp
- Pull Placement Methods

Sunny Medical Equipment Limited

Room1208,Jingying Building,No.316,Huanshi Middle Road,Xiaobei,Yuexiu,Guangzhou,China
 Cel&Whatsapp: 0086-18680207678 Email: sale008@sunnymedeqp.com

| Product No. | | Size | O.D.(mm) | Length(mm) | Remark |
|--------------|---------------|------|----------|------------|----------------------------|
| Economic Kit | Luxury kit | | | | |
| SY-L147-16SE | SY-L147A-16SL | 16FR | 5.3 | 600 | Length could be customized |
| SY-L147-18SE | SY-L147A-18SL | 18FR | 6.0 | | |
| SY-L147-20SE | SY-L147A-20SL | 20FR | 6.7 | | |
| SY-L147-22SE | SY-L147A-22SL | 22FR | 7.3 | | |
| SY-L147-24SE | SY-L147A-24SL | 24FR | 8.0 | | |

PEG

PEG Kit (Luxury Pack)

Traction Removable PULL PEG Tube × 1
 5 Fr Retrieval Snare (1.8mm × 160cm) × 1
 Feeding Adapter × 1
 Feeding Adapter with Luer lock × 1
 Retention Ring × 1
 25Ga Needle × 1
 19/21Ga Needle × 1
 Stainless Steel Scissors × 1
 Stainless Steel Hemostat × 1
 Tubing Clamp × 1
 Water Soluble Lubricant × 1
 10*10 Gauze Sponges × 4
 5*5 Slit Gauze Sponges × 2
 Placement Wire (0.36/0.6mm × 260cm) × 1
 10ml Syringe × 1
 Disposable Scalpel with #11 Blade × 1
 Introducer Needle × 1
 Fenestrated Drape × 1
 Directions for Use × 1

PEG

PEG Kit (Economic Pack)

Traction Removable PULL PEG Tube × 1
 Feeding Adapter × 1
 Feeding Adapter with Luer lock × 1
 Retention Ring × 1
 19/21Ga Needle × 1
 25Ga Needle × 1
 Placement Wire (0.36/0.6mm × 260cm) × 1
 10ml Syringe × 1
 Disposable Scalpel with #11 Blade × 1
 Introducer Needle × 1
 Tubing Clamp × 1
 10*10 Gauze Sponges × 4
 Fenestrated Drape × 1
 Directions for Use × 1

