



From Indigenous Ideas to Newtech Solutions



# ClearEmbo

*Embolectomy Catheter*

## Newtech Medical Device

We develop & manufacture, market and support cutting edge vascular access devices and accessories to meet the clinical needs of the medical industry, particularly in the field of interventional cardiology and dialysis. Our company's engineering and applications expertise provides superior products whose progressive designs accommodate advances in medical device and whose quality anticipates the requirements of our professional clients and the patients they serve.

Currently one of the largest manufacturer & importer in the field of Cardiology, Nephrology & Surgical Disposables in INDIA.

## Embolectomy Catheter

Embolectomy catheters are sterile medical devices dedicated for removal of fresh arterial embolism. Catheter is non-toxic, pyrogene free, visible under fluoroscopy.

Single lumen catheters use hand-tied balloons that are glued in three areas and pre-stretched for flexibility. These catheters are strong and able to withstand a substantial force before the balloon disengages from the catheter. Balloons are well-centered and resistant to rupture.

## Product Feature

### Resistant to rupture

Balloon of our Embolectomy catheters are manufactured in this manner that it is resistant to rupture during process.

### Strength

Constructed with durable material which help to withstand a substantial force before the balloon disengages.



### Colour coding and pacing

Embolectomy Catheters are colour coded to visually segregate various French Sizes. These are packed in strong plastic shell type casing.



### Predictable inflation

The balloons inflate at a consistent rate as the volume of the balloon increases



### With Stylet

Newtech Embolectomy Catheters comes with stainless stylet for easier insertion of the catheter.

## Specifications

REF	French	Length	Diameter of Inflated Balloon		Colour Code
			φ mm	0.9% NaCl(ml)	
NT/EC/SL240	02 Fr	40 cm	4	0.15	Transparent
NT/EC/SL260	02 Fr	60 cm*	4	0.15	Transparent
NT/EC/SL380	03 Fr	80 cm*	6	0.2	Green
NT/EC/SL480	04 Fr	80 cm*	8	0.7	Red
NT/EC/SL580	05 Fr	80 cm*	10	1.5	White
NT/EC/SL680	06 Fr	80 cm*	12	2.0	Blue
NT/EC/SL780	07 Fr	80 cm*	14	2.5	Yellow
NT/EC/SL880	08 Fr	80 cm*	16	3.5	Brown
NT/EC/SL980	09 Fr	80 cm	18	4.0	Purple
NT/EC/SL1080	10 Fr	80 cm	20	4.5	Grey

No Star - No Sampling, Check Availability before confirming order of Non Running products.

\* No Sampling Available, Check Availability of Running Products.

\*\* Sampling Available, No need to check availability of highly moving products, Easily Available

## Newtech Medical Devices,

14/5 3rd Floor , Mathura Road , Near Mewla Maharajpur  
Metro Station, Faridabad, Haryana -121003, INDIA

Tel +91 129 225 9961

+91 129 225 9962

Mail ID info@ntmdevices.com

Website www.ntmdevices.com

Our Distributor :