

# Mediblu

INTEGRATED  
**MEDICAL**  
SOLUTIONS

BRO-Mediblu-Rev D-October 2021

3016 NW 82 Ave Doral FL, 33122, United States | +1 302 451.92.34

[www.mediblumed.com](http://www.mediblumed.com)

# PATIENT MONITORING

## MM3

VITAL SIGNS MONITOR



- » SpO2+NIBP+ TEMP
- » 5.6" high resolution color TFT-LCD screen Lightweight & portable design
- » User-friendly interface
- » Flexible configurations to meet different clinical needs Auto/Manual/Continuous/Average BP mode Real-time measurements & trend display
- » Built-in rechargeable Lithium-ion battery for 8 hours of continuous work
- » USB flash disk support

## MM1

VITAL SIGNS MONITOR



- » 7" High resolution (1024\*600) colour TFT display
- » Suitable for Clinics, GP Surgeries, Hospitals, EMS, Patient-Transfer, Ambulance.
- » Suitable for adults, pediatric, neonates.
- » Data memory storage of up to 360 hours.
- » High-capacity lithium battery for continuous patient monitoring.
- » Optional ETCO2, Touchscreen, Nellcor SpO2
- » Portable Design
- » Bluetooth ear thermometer
- » USB output

## MM12

MULTI-PARAMETER PATIENT MONITOR



- » 12.1" color LCD multichannel waveform display
- » SpO2, ECG, NIBP, TEMP, RR. » Optional: IBP, ETCO2, C.O
- » Rechargeable high-energy built-in battery.
- » Easy operation with user friendly menu structure
- » Defibrillator protection
- » Central Monitoring
- » Thermal printer and touch screen (optional).
- » 360 hours data storage.
- » Drug dose calculation function.
- » Suitable for adult, pediatrics, neonates

## Blu12

MULTI-PARAMETER PATIENT MONITOR



- » 12.1" Touch Screen (optional) colour LCD screen with 8-waveform display.
- » Flat Display with adjustable lightness and volume to improve simultaneous display of multi-channel real time waveform & trend diagram.
- » Battery: Rechargeable high-energy built-in battery.
- » Operating menu with multi-language interface selection.
- » User friendly menu structure design, touch screen and rotary dial.
- » Special protections against the ESU interference & defibrillation, no need to disconnect the monitor from the patient during defibrillation.
- » Three-application modes: monitoring, diagnosis, operating.
- » 72 hours data storage.
- » Intelligent audio-visual comprehensive alarm.

## Blu12+

MULTI-PARAMETER PATIENT MONITOR



- » 12.1" Touch Colour LCD Screen.
- » Multi-channel waveform display.
- » Portable operation, multi-pronged.
- » Full-screen touch, fingertip operation.
- » Fast button backlit design for easy night operation.
- » Ergonomic knob operation, classic and flexible.
- » Plugin & use, powerful 2 module slots, easy to upgrade.
- » Flexible combination, resource sharing.
- » Global gold standard Masimo SpO2, ensures low perfusion measurement.
- » Sun tech NIBP blood pressure, accurate blood pressure measurement.

## Blu15+

MULTI-PARAMETER PATIENT MONITOR



- » 15" Touch Colour LCD Screen.
- » Multi-channel waveform display.
- » Portable operation, multi-pronged.
- » Full-screen touch, fingertip operation.
- » Fast button backlit design for easy night operation.
- » Ergonomic knob operation, classic and flexible.
- » Plugin & use, powerful 3 module slots, easy to upgrade.
- » Flexible combination, resource sharing.
- » Global gold standard Masimo SpO2, ensures low perfusion measurement.
- » Sun tech NIBP blood pressure, accurate blood pressure measurement.

**FD1**

ULTRASONIC POCKET DOPPLER



- » Compact design with easy one-hand operation
- » Hi-Fi speaker
- » Interchangeable Waterproof Probe
- » 2 MHz probe for later pregnancy and larger patients
- » 3MHz probe fetal heartbeat as early as 9 weeks
- » 4/5/8 MHz probe for blood flow
- » Detection

**FM02**

FETAL MATERNAL MONITOR



- » Compact, portable and lightweight
- » Long life Lithium-ion battery for outpatient service
- » 5.6" foldable screen
- » Large numeric and waveform display
- » Simple workflow with silicone buttons for operations
- » Signals Overlap Verification to differentiate twins FHR
- » FHR signal quality indicator helps optimize the probe position
- » MFM-CNS network for remote monitoring
- » Insight software for data management on PC

**FM10**

FETAL MONITOR



- » 10.4 inch color TFT LCD screen
- » Calculation of auto/manual movement
- » Auto/manual patient ID generation, easy user management
- » Three monitoring modes: fetal, fetal and maternal, maternal
- » Real-time printing and timing printing
- » Built-in high capacity rechargeable battery
- » Wired, wireless and hybrid central monitoring solution.
- » FHR + TOCO + Fetal movement + maternal ECG + NIBP + Temperature
- » + pulse rate + SpO2 + respiration + battery +110mm printer

**DVC300**

VIDEO COLPOSCOPE



- » HD Option
- » Captures details with higher accuracy and with extended magnification (46X)
- » HD camera
- » Innovative 2,5" screen display.
- » Complete and versatile diagnostic information critical in clinical decision.
- » A variety of image modes, True tissue colour, electronic green filter
- » Unique light source design
- » Auto-focus control with AUTO adjustment of white balance
- » 360° adjustable full-damping camera holder

**HA3P**


ECG HOLTER SYSTEM



- » Lightweight design: only 50 gr
- » High level of waterproof
- » 3 or 12-channel digital recording
- » Real-time waveforms display
- » Up to 8 days recording with one AAA battery
- » Complete Software for Data Processing
- » GDT/XML protocol support

**S12/S12W**

PC ECG / Stress Test



- » Ability to work with both PC & Tablet
- » Wireless Solution
- » 12-channel Resting ECG
- » Stress test
- » Comprehensive Analysis
- » Versatile Report PDF, FDA-XML, SCP, DICOM, JPG, BMP and DOC
- » Format reports transferable via Email
- » Compatible with GDT Systems CG, Frequency spectrum cardiogram
- » Compatible with reputable treadmills and bicycle

# CARDIOLOGY

## ME3

MULTI CHANNEL  
ELECTROCARDIOGRAPH



- » 3-channels ECG
- » 5.1 inch high-resolution LCD
- » Durable and comfortable silicon keys
- » Rechargeable lithium battery for 2.5 hours operation
- » Peripheral ports: SD card, USB Up to 20 files internal storage
- » Unlimited files with external SD card
- » USB port for data transmission to PC

## ME6P

MULTI CHANNEL  
ELECTROCARDIOGRAPH



- » 6-channels ECG
- » 8 inch TFT LCD touching screen
- » Durable and comfortable silicon keys and knob
- » Rechargeable lithium battery for 2.5 hours operation
- » Peripheral ports: SD card, USB, LAN
- » File Management Software
- » Up to 200 files internal storage

## ME12P

MULTI CHANNEL  
ELECTROCARDIOGRAPH



- » 12-channels ECG
- » 8.9 inch TFT LCD touch screen
- » 60 degree viewing angle adjustable
- » Full alphanumeric waterproof silicon keyboard
- » Abundant Peripheral interfaces: VGA, LAN, USB, SD CARD
- » ECG data management software
- » High resolution Printer A4 paper size
- » Storage: 1500 ECG Files
- » Rechargeable Li-ion battery, 14.8V-4400mAh

## M01

PULSE OXIMETER



- » Real-time SpO2 waveform.
- » SpO2 & PR value and intensity.
- » Fast, accurate, stable & strong anti-interference performance.
- » At least 48 hours parameter & trend data storage, up to 999 measuring objects.
- » PC connection to download & analyze/print data.
- » Suitable for adult, pediatric & neonatal use.
- » Real time clock.
- » Automatically recharges battery upon circuit connection.
- » 10-12 hours battery back-up, Real capacity indication.
- » Visual and acoustic Abnormal Parameter Alarm.

**Mediblu**

INTEGRATED  
**MEDICAL SOLUTIONS**

3016 NW 82 Ave Doral FL, 33122, United States

+1 305 406 39.31

+1 305 406 21.44



[www.mediblumed.com](http://www.mediblumed.com)